



TERUMO

CAPIOX[®] BUBBLE TRAP

PIEGE A BULLES / BLASENFALLE / ATRAPA BURBUJAS
TRAPPOLA PER BOLLE / BUBBLE TRAP

CE 0197



KEEP AWAY
FROM SUNLIGHT

STER/ESTER:

LOT:

STERILE EO

Sterilized by ethylene oxide.



Read the instructions for use. For single use only.

CAUTION: Federal law(U.S.A.) restricts this device to sale by or on the order of a physician.



TERUMO CORPORATION TOKYO 151-0072, JAPAN MADE IN JAPAN

TERUMO CARDIOVASCULAR SYSTEMS ANN ARBOR, MICHIGAN 48103-9300, U. S. A.

EC REP TERUMO EUROPE N. V. 3001, LEUVEN, BELGIUM


1

UNIT
UNITE
STÜCK

UNIDAD
UNITA
STUK



KEEP AWAY
FROM RAIN



TEMPERATURE
LIMITS



(01)04987350705730(17)160100(10)130205

®:Registered Trademark.

CAPIOX® Bubble Traps

	BT05	BT15
Blood port	1/4" (6.4 mm)	3/8" (9.5 mm)
Vent port	Female luer	Female luer
Filter material	Polyester	Polyester
Filter pore size	170 microns	170 microns
Housing	Polycarbonate	Polycarbonate
Prime volume	50 mL	150 mL