

Carefully Crafted



SurGuard[®] 3

Safety Hypodermic Needle

Safety made easy


- Supports all 3 ways of activation
- Audible, visible and tangible confirmation²
- User comfort – leading to patient's confidence²
- EN ISO 23908 compliant

Drug Delivery Devices

Features and benefits

- Hinged safety sheath – stable, reliable, easy to adopt²
- Sharper needles for more comfortable injections¹
- Sheath locks to both the needle and the hub to provide added security and visible, audible and tactile feedback²
- Finger & thumb grips designed to suit all hand sizes
- Can be used on luer slip tip and luer lock syringes
- 3 ways of shield activation to suit individual clinician needs and preferences: hard surface or finger or thumb

Sheath activation



¹ Comparative Hypodermic Needle Sharpness Among Major Brands: Results of an Independent Study 2009

² Data on file

Please contact us for more details.

Refer to the product labels and instructions for use for warnings, cautions and precautions.

SurGuard®3 Safety Hypodermic Needle

References	Hub Color	Diameter (OD)		Length	
		Gauge	mm	mm	Inch
SG3+1813	Pink	18	1.2	13	1/2
SG3+1825	Pink	18	1.2	25	1
SG3+1838	Pink	18	1.2	38	1 1/2
SG3+1916	Ivory	19	1.1	16	5/8
SG3+1938	Ivory	19	1.1	38	1 1/2
SG3+2025	Yellow	20	0.9	25	1
SG3+2038	Yellow	20	0.9	38	1 1/2
SG3+2051	Yellow	20	0.9	51	2
SG3+2125	Green	21	0.8	25	1
SG3+2138	Green	21	0.8	38	1 1/2
SG3+2151	Green	21	0.7	51	2
SG3+2213	Black	22	0.7	13	1/2
SG3+2225	Black	22	0.7	25	1
SG3+2232	Black	22	0.7	32	1 1/4
SG3+2238	Black	22	0.7	38	1
SG3+2251	Black	22	0.7	51	2
SG3+2325	Blue	23	0.6	25	1
SG3+2338	Blue	23	0.6	38	1 1/2
SG3+2516	Orange	25	0.5	16	5/8
SG3+2525	Orange	25	0.5	25	1
SG3+2609	Brown	26	0.45	9	3/8
SG3+2713	Grey	27	0.4	13	1/2
SG3+3013	Yellow	30	0.30	13	1/2

Please note that some products may not be available in all regions, check with your local TERUMO Pharmaceutical Solutions contact for details.

sterilized using irradiation **CE 0197**
The device has no components made of natural rubber latex or PVC.