



## Radifocus® Angiographic Catheter



[Enlarge image](#)



Radifocus® Angiographic Catheter is a 5 Fr. angiographic catheter intended for use in angiographic procedures.

It delivers radiopaque media and therapeutic agents to selected sites in the vascular system.

It is also used to lead a guide wire or a catheter into the target site.

### Product Characteristics

- Flow, smoothness and superior safety thanks to the polyamide surface, the round tip and one-piece body structure
- Soft and non traumatic catheter with excellent torque control due to its internal metallic mesh braid
- Smooth radiopaque surface to prevent thrombus formation
- High pressure resistance (up to 1000 psi / 6.895 kPa)
- Various tip shapes available

### General Specifications

Outer diameter	5 Fr / 1.70 mm
Inner diameter	0.043" / 1.10 mm
Maximum guide wire size	0.038" / 0.97 mm
Pressure limit	1000 psi / 6,895 kPa
Catheter material	Polyamide
Mesh braid material	Stainless steel

### Item Specifications

Catheter length	Distal curve	Side hole	Item reference
65 cm	<a href="#">Pigtail</a>	6	RF-FG1500GM
80 cm	<a href="#">Pigtail</a>	6	RF-FG15008M
110 cm	<a href="#">Pigtail</a>	6	RF*FG15011M
65 cm	<a href="#">Straight</a>	6	RF-FK1500GM
80 cm	<a href="#">Straight</a>	6	RF-FK15008M
80 cm	<a href="#">Straight 8 cm tapered</a>	0	RF*DI15008M
80 cm	<a href="#">Renal</a>	2	RF-DA15008M
80 cm	<a href="#">Renal</a>	0	RF-DA25008M
80 cm	<a href="#">Renal Left</a>	0	RF-DAC5008M
65 cm	<a href="#">Cobra Small</a>	2	RF-DB1500GM
65 cm	<a href="#">Cobra Small</a>	0	RF-DB4500GM
80 cm	<a href="#">Cobra Small</a>	0	RF-DB45008M
65 cm	<a href="#">Cobra Middle</a>	2	RF-DB2500GM
65 cm	<a href="#">Cobra Middle</a>	0	RF-DB5500GM
80 cm	<a href="#">Cobra Middle</a>	0	RF-DB55008M
65 cm	<a href="#">Cobra Large</a>	2	RF-DB3500GM
65 cm	<a href="#">Cobra Large</a>	0	RF-DB6500GM
80 cm	<a href="#">Cobra Large</a>	0	RF-DB65008M
65 cm	<a href="#">Shepherd Hook 0.8 cm</a>	0	RF-DC3500GM
65 cm	<a href="#">Shepherd Hook 1.0 cm</a>	0	RF-DC4500GM
65 cm	<a href="#">J curve Small</a>	0	RF-DD4500GM
65 cm	<a href="#">J curve Middle</a>	0	RF-DD5500GM
65 cm	<a href="#">J curve Large</a>	0	RF-DD6500GM
100 cm	<a href="#">Simmons 1</a>	0	RF-EA15010M
100 cm	<a href="#">Simmons 2</a>	0	RF-EA25010M
100 cm	<a href="#">Simmons 3</a>	0	RF-EA35010M
100 cm	<a href="#">Headhunter 1</a>	0	RF-EB15010M
100 cm	<a href="#">Headhunter 3</a>	0	RF-EB25010M
100 cm	<a href="#">Bentson-Hanafee-Wilson 2</a>	0	RF-EE25010M
100 cm	<a href="#">MANI</a>	0	RF-EG15010M
65 cm	<a href="#">Vertebral</a>	0	RF-EH1500GM
100 cm	<a href="#">Vertebral</a>	0	RF*EH15010M
100 cm	<a href="#">Hinck 1</a>	0	RF-EJ15010M

Please quote above item reference codes when placing an order.

Other code numbers are available on special demand. For any further information, please contact your local [Terumo representative](#).