

# LOOKS SHARP. FEELS SOFT.

Insulin delivery with  
Japanese care



 **TERUMO**

**FINEGLIDE™**  
SURPRISINGLY SOFT



## INSULIN PEN NEEDLE SOLUTION FOR SELF-INJECTION.



### 3-BEVEL SUPER SHARP

With High-Grade Silicone surface, the injection is more gentle, comfortable and pain-free.



### MAXIMUM FLOW RATE

Thin-wall technology enhances the medication flow providing a comfortable experience.



### ONE FOR ALL "UNIVERSAL FIT"

Designed to fit all insulin pens and GLP-1 RECEPTOR AGONIST PENS\* standard ISO 11608.



**100 YEARS OF TRUST**

We are committed to contributing to society through healthcare for 100+ years. Headquartered in Tokyo and operating globally with a presence in more than 160 countries and regions through our strong network of 26000+ associates, we help to provide innovative medical care solutions.



**Thin-Wall technology, reduced injection force providing a comfortable experience.**

**PRODUCT LINEUP**

**NEEDLE LENGTH**

Available in 29G to 32G with varying needle lengths from 4mm upto 12mm



**31G**  
4mm



**32G**  
4mm



**32G**  
5mm







**31G**  
6mm



**32G**  
6mm

Average injection force for **FINEGLIDE™**: 31% better than the needles available in market



Packing	100 needles in a box   20 pouches of 5 needles each	
Sterilization	Ethylene Oxide	
Certificate	ISO 11608-2    ISO 13485	  
Material	Needle: Stainless steel (Medical Grade) Needle hub: Polypropylene	Needle cap: Polyethylene Primary container: Polypropylene

**FINEGLIDE™**  
SURPRISINGLY SOFT

TERUMO INDIA PRIVATE LIMITED.

16TH FLOOR, TOWER B, UNITECH CYBER PARK, SECTOR-39, GURUGRAM- 122001, HARYANA, INDIA.

TEL: +91.124.4718700, FAX: +91.124.4718718 CIN:U33110HR2013FTC049841